



# Health and Medical Supply




# OUR PRODUCTS


HealthandMedicalSupply.com LLC is a Florida-based small business dedicated to helping get medical products to those who need them with the highest professional and ethical standards.

## MASKS


## N95-MASKS




**Alpha Pro Tech**  
**N-95 Particulate Respirator**  
PFL® (Positive Facial Lock®) Design  
with Magic Arch® and Twist Seal®




**Donning Instructions**




**Step 1.** Grasp top of mask with one hand. Slide top edge of mask over forehead and over nose.




**Step 2.** Open mask by pulling lower edge and use other hand to pull mask to fully open. If in NIOSH (N95) that mask is fully open.




**Step 3.** Place over face in pocket position. Bend back over back of head, position, behind head behind neck behind ears and top band over crown of head.




**Step 4.** Cover nose and mouth to ensure full contact, making it the Magic Arch® to fully extend along your nose and mouth.




**Step 5.** Adjust mask to ensure proper fit. Adjust the mask on the forehead above the eyes, and ensure properly fit that mask in NIOSH approved form.



**Step 6.** Check any necessary sealant together under chin and ear over head.



**Step 7.** Final Check: To check fit, place hand against the top edge of the respirator and inhale. If air escapes around your nose, top of the respirator or under the chin, repeat steps 1 and 2.



**Step 8.** PFL. Look closely at face.

**Warnings:**  
1. This respirator should be used only for the purposes for which it is designed. Change filters and/or cartridges when they become clogged or saturated. Do not use this respirator if the filter or cartridge is damaged or if the seal is broken.  
2. This respirator is not intended for use in environments where there are airborne infectious agents or other hazardous substances.  
3. This respirator is not intended for use in environments where there are airborne radioactive materials or other hazardous substances.  
4. This respirator is not intended for use in environments where there are airborne chemical agents or other hazardous substances.  
5. This respirator is not intended for use in environments where there are airborne biological agents or other hazardous substances.

**Notes:**  
A. Before use, the respirator and its accessories should be inspected to ensure that they are in proper condition and ready for use.  
B. The respirator should be stored in a clean, dry, and well-ventilated area.  
C. The respirator should be stored in a clean, dry, and well-ventilated area.  
D. The respirator should be stored in a clean, dry, and well-ventilated area.  
E. The respirator should be stored in a clean, dry, and well-ventilated area.

**AlphaProTech**  
LLC

Customer Service - Phone: 800.749.1363 • Fax: 905.479.8783 • sales@alphaprotech.com



## **ALPHA PROTECH CUSTOM FIT NIOSH MASK MODEL 695**

**NIOSH Approved**, our N-95 delivers the appropriate amount of filtration efficiency and when combined with a trio of patented and unique design features, also delivers the user a comfortable, custom fitting medical facemask. The Positive Facial Lock® (PFL®) N-95 Particulate Respirator meets the World Health Organization (WHO) and the Centers for Disease Control (CDC) recommended protection levels (where applicable) for many airborne contaminants.

Terms are prepaid unless otherwise indicated or determined on a case by case basis. Financing is available for most products.

[Inquire for distributor and government pricing](#)

**SKU:** HMS-APT-312

**Price:** \$2.95

**Categories:** [Masks](#), [N95-Masks](#), [Our Products](#)





# **Health and Medical Supply**

