



Health and Medical Supply



COVID TESTS



OUR PRODUCT

10 MINUTES
DETECTION TIME

• Just one finger prick of blood



ONLY €X

- Customs
 - Insurance
 - Shipping
 - Freight Excluded
- Only Port Delivery if Turkey preferred

2

HOW IT WORKS?

CURRENT METHOD




WHO urges world to "test, test, test" for COVID-19

3

HOW IT WORKS?

PRE-SCREENING



BENEFITS of pre-screening

- Exponentially **faster**
- **No Queues, No Wait** time for results
- Exponentially more **cost effective**
- **Early Detection & Elimination**
- Effort less, and **accessible** to everyone nationwide, instantly
- Viable **precautionary** step

HOW IT WORKS?


Detection Period 10 MINUTES

- **Easy operation** without requirement of any Doctor or Professional Nurse
- No special equipment storage and transportation conditions required
- Works with whole blood, serum, and plasma
- **Tests for 2 antibodies** IgM and IgG simultaneously
- **Instant field screening**



98%+ACCURACY

HOW IT WORKS?



| | Fluorescence PCR | CUA | Colloidal Gold Method (our method) |
|---------------------------------|--|--------------------------------|------------------------------------|
| Detection substance | Nucleic acid | Antibody | Antibody |
| Type of sample | Nasopharyngeal swabs, sputum, saliva, swabs, nasal | Serum/ Plasma | Serum/ Plasma / Whole blood |
| Time to get result | 2 hrs | 30 min | Fastest 10 minutes |
| Instrument needed or not | Yes | Yes (3000, 1 1000) | Not needed |
| Laboratory requirement | High | Relatively High | Average |
| Product usage | Confirming diagnosis | Nucleic acid negative/positive | Nucleic acid negative/positive |

Colloidal Gold Method: The 2019-nCoV IgG/IgM Rapid Test Device using this method is a rapid chromogenic immunoassay for the qualitative detection of IgG & IgM antibody of Coronavirus in human whole blood, serum, or plasma.



RT-PCR Vs SARS-CoV-2 IgM Ab Rapid Test Comparative Test Summary

In this trial, 1300 clinical samples were selected, there were 300 positive samples and 1000 negative samples.

The SARS-CoV-2 IgM Ab rapid test and the SARS-CoV-2 RT-PCR test were detected simultaneously and the Sensitivity (positive coincidence rate), Specificity (negative coincidence rate), and Accuracy (total coincidence rate) were calculated:

| Antibody | Sensitivity | Specificity | Accuracy |
|----------|-------------|-------------|----------|
| IgG | 93% | 97.5% | 96.5% |
| IgM | 92% | 96% | 95.8% |

COVID 19 TOTAL ANTIBODY TEST (EUROPEAN FDA EUA SUBMITTED)

Covid- 19 total antibody, IgM/IgG Combo Rapid Test. For the simultaneous detection of IgM and IgG antibodies to SARS-CoV-2 virus in human whole blood, serum or plasma.

Origin: US, **Wait Time:** 2-3 weeks, **Minimum Order Quantity:** 10K


Terms are prepaid unless otherwise indicated or determined on a case by case basis. Financing is available for most products.

[Inquire for distributor and government pricing](#)

SKU: C19-RAP-TST-MX

Price: \$10.90

Categories: [Covid Tests](#), [Our Products](#)



QSARS-COV-2 IGG/IGM RAPID TEST

qSARS-CoV-2 IgG/IgM Cassette Rapid Test is a product for in-vitro analysis of whole blood, serum, and plasma, designed to deliver qualitative results for a panel of tests.

Origin: US, **Wait Time:** 1-2 Weeks, **Minimum Order Quantity:** 10K

Terms are prepaid unless otherwise indicated or determined on a case by case basis. Financing is available for most products.

[Inquire for distributor and government pricing](#)

SKU: C19-RAP-TST-CI

Price: \$27.30

Category: [Covid Tests](#)



COVID-19 IGG/IGM RAPID TEST CASSETTE (WHOLE BLOOD/SERUM/PLASMA)

The COVID-19 IgG/IgM (Whole Blood/Serum/Plasma) Rapid Test Device utilizes lateral flow technology that is used for the qualitative, differential detection of both anti- SARS-CoV-2 IgM and IgG antibodies. In general, antibodies can be detected 1-3 weeks after infection. This test is intended to screen patients for COVID-19. Combining RNA and Antibody tests can significantly raise the sensitivity for detecting COVID-19 in infected individuals.

Origin: US, **Wait Time:** 1-2 Weeks, **Minimum Order Quantity:** 10K

Terms are prepaid unless otherwise indicated or determined on a case by case basis. Financing is available for most products.

[Inquire for distributor and government pricing](#)

SKU: C19-RAP-TST-HG

Price: \$26.00

Category: [Covid Tests](#)



Health and Medical Supply

